

产品名称: Serpin A11 Rabbit Polyclonal Antibody
产品货号: APRab17763



产品概述 (Summary)

产品名称 (Production Name)	Serpin A11 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SERPINA11
别名 (Alternative Names)	SERPINA11; Serpin A11
基因 ID (Gene ID)	256394.0
蛋白 ID (SwissProt ID)	Q86U17.The antiserum was produced against synthesized peptide derived from human SERPINA11. AA range:191-240

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	47kDa

研究背景 (Background)

产品名称: Serpin A11 Rabbit Polyclonal Antibody
产品货号: AP Rab17763



similarity:Belongs to the serpin family.,similarity:Belongs to the serpin family.,

研究领域 (Research Area)

Immunology; Adaptive Immunity; B Cells; Non-CD; Cancer; Oncoproteins/suppressors; Oncoproteins; Cell Biology; Proteolysis / Ubiquitin; Protease inhibitors; Serine protease inhibitors; SERPINs

图片 (Image Data)



Western blot analysis of lysates from HT-29 cells, using SERPINA11 Antibody. The lane on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .