

产品名称: SFRS15 Rabbit Polyclonal Antibody  
产品货号: APRab17799



## 产品概述 (Summary)

|                        |                                   |
|------------------------|-----------------------------------|
| 产品名称 (Production Name) | SFRS15 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody        |
| 宿主 (Host)              | Rabbit                            |
| 应用 (Application)       | WB,IHC,ICC/IF,ELISA               |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                   |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                            |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | SCAF4                                                                                                                                      |
| 别名 (Alternative Names) | SCAF4; KIAA1172; SFRS15; Splicing factor; arginine/serine-rich 15; CTD-binding SR-like protein RA4; SR-related and CTD-associated factor 4 |
| 基因 ID (Gene ID)        | 57466.0                                                                                                                                    |
| 蛋白 ID (SwissProt ID)   | O95104.The antiserum was produced against synthesized peptide derived from human SFRS15. AA range:601-650                                  |

## 产品应用 (Application)

|                          |                                                                         |
|--------------------------|-------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000 |
| 蛋白分子量 (Molecular Weight) | 140kDa                                                                  |

产品名称: SFRS15 Rabbit Polyclonal Antibody  
产品货号: APRab17799




## 研究背景 (Background)


This gene likely encodes a member of the arginine/serine-rich splicing factor family. A similar protein in Rat appears to bind the large subunit of RNA polymerase II and provide a link between transcription and pre-mRNA splicing. Alternatively spliced transcript variants have been described. [provided by RefSeq, Feb 2009],function:May act to physically and functionally link transcription and pre-mRNA processing.,similarity:Belongs to the splicing factor SR family.,similarity:Contains 1 CID domain.,similarity:Contains 1 RRM (RNA recognition motif) domain.,subunit:Interacts with the repetitive C-terminal domain (CTD) of RNA polymerase II.,

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of lysates from COLO cells, using SFRS15 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of COLO cells using SFRS15 Polyclonal Antibody diluted at 1: 2000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA) .

## 注意事项 (Note)

For research use only .

产品名称: SFRS15 Rabbit Polyclonal Antibody  
产品货号: APRab17799

---

