

产品名称: SFT2D2 Rabbit Polyclonal Antibody
产品货号: APRab17802



产品概述 (Summary)

产品名称 (Production Name)	SFT2D2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SFT2D2
别名 (Alternative Names)	SFT2D2; Vesicle transport protein SFT2B; SFT2 domain-containing protein 2
基因 ID (Gene ID)	375035.0
蛋白 ID (SwissProt ID)	O95562.The antiserum was produced against synthesized peptide derived from human SFT2B. AA range:61-110

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	18kDa

研究背景 (Background)


产品名称: SFT2D2 Rabbit Polyclonal Antibody
产品货号: APRab17802



function: May be involved in fusion of retrograde transport vesicles derived from an endocytic compartment with the Golgi complex., similarity: Belongs to the SFT2 family., function: May be involved in fusion of retrograde transport vesicles derived from an endocytic compartment with the Golgi complex., similarity: Belongs to the SFT2 family.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from LOVO cells, using SFT2B Antibody. The lane on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .