

产品名称: SH2D4A Rabbit Polyclonal Antibody
产品货号: APRab17836



产品概述 (Summary)

产品名称 (Production Name)	SH2D4A Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SH2D4A
别名 (Alternative Names)	SH2D4A; SH2A; SH2 domain-containing protein 4A; Protein SH(2)A
基因 ID (Gene ID)	63898.0
蛋白 ID (SwissProt ID)	Q9H788. The antiserum was produced against synthesized peptide derived from human SH2D4A. AA range: 180-220

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)

产品名称: SH2D4A Rabbit Polyclonal Antibody
产品货号: AP Rab17836




sequence caution:Translated as stop.,sequence caution:Translation N-terminally extended.,similarity:Contains 1 SH2 domain.,tissue specificity:Ubiquitously expressed. Aberrantly expressed in some cancers.,sequence caution:Translated as stop.,sequence caution:Translation N-terminally extended.,similarity:Contains 1 SH2 domain.,tissue specificity:Ubiquitously expressed. Aberrantly expressed in some cancers.,

研究领域 (Research Area)

Tags & Cell Markers; Cell Type Markers; Tumor Associated; Cytoskeleton / ECM; Cell Adhesion; Cell Adhesion Molecules; Endothelial; Signal Transduction; Protein Phosphorylation; Tyrosine Kinases

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200 (4° overnight) . 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200 (room temperature, 30min) .

注意事项 (Note)

For research use only .