

产品名称: SH2D5 Rabbit Polyclonal Antibody
产品货号: APRab17837



产品概述 (Summary)

产品名称 (Production Name)	SH2D5 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SH2D5
别名 (Alternative Names)	SH2D5; SH2 domain-containing protein 5
基因 ID (Gene ID)	400745.0
蛋白 ID (SwissProt ID)	Q6ZV89.The antiserum was produced against synthesized peptide derived from human SH2D5. AA range:281-330

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	35kDa

研究背景 (Background)


产品名称: SH2D5 Rabbit Polyclonal Antibody
产品货号: APRab17837




similarity:Contains 1 SH2 domain.,similarity:Contains 1 SH2 domain.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of SH2D5 Antibody. The lane on the right is blocked with the SH2D5 peptide.



Western Blot analysis of various cells using SH2D5 Polyclonal Antibody

注意事项 (Note)

For research use only .