

产品名称: SH3BGR Rabbit Polyclonal Antibody
产品货号: APRab17839



产品概述 (Summary)

产品名称 (Production Name)	SH3BGR Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SH3BGR
别名 (Alternative Names)	SH3BGR; SH3 domain-binding glutamic acid-rich protein; SH3BGR protein; 21-glutamic acid-rich protein; 21-GARP
基因 ID (Gene ID)	6450.0
蛋白 ID (SwissProt ID)	P55822.The antiserum was produced against synthesized peptide derived from human SH3BGR. AA range:81-130

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	26kDa

产品名称: SH3BGR Rabbit Polyclonal Antibody
产品货号: AP Rab17839




研究背景 (Background)

caution:It is uncertain whether Met-1 or Met-64 is the initiator.,similarity:Belongs to the SH3BGR family.,tissue specificity:Expressed in heart and skeletal muscle.,caution:It is uncertain whether Met-1 or Met-64 is the initiator.,similarity:Belongs to the SH3BGR family.,tissue specificity:Expressed in heart and skeletal muscle.,

研究领域 (Research Area)

Signal Transduction; Adapters; Cytoplasmic; Tags & Cell Markers; Cell Type Markers; Other Cell Types

图片 (Image Data)



Western blot analysis of lysates from Jurkat cells, using SH3BGR Antibody. The lane on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .