

产品名称: SLC10A7 Rabbit Polyclonal Antibody
产品货号: APRab17951



产品概述 (Summary)

产品名称 (Production Name)	SLC10A7 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SLC10A7
别名 (Alternative Names)	SLC10A7; C4orf13; P7; PSEC0051; Sodium/bile acid cotransporter 7; Na(+)/bile acid cotransporter 7; Solute carrier family 10 member 7
基因 ID (Gene ID)	84068.0
蛋白 ID (SwissProt ID)	Q0GE19.The antiserum was produced against synthesized peptide derived from human SLC10A7. AA range:168-217

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	39kDa

产品名称: SLC10A7 Rabbit Polyclonal Antibody
产品货号: APRab17951




研究背景 (Background)


similarity:Belongs to the sodium:bile acid symporter family.,similarity:Belongs to the sodium:bile acid symporter family.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from mouse liver, using SLC10A7 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from 293 cells using SLC10A7 antibody.

注意事项 (Note)

For research use only .