

产品名称: SLC17A2 Rabbit Polyclonal Antibody
产品货号: AP Rab17953



产品概述 (Summary)

产品名称 (Production Name)	SLC17A2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SLC17A2 SLC17A2; NPT3; Sodium-dependent phosphate transport protein 3; Na(+)/PI cotransporter 3; Sodium/phosphate cotransporter 3; Solute carrier family 17 member 2
别名 (Alternative Names)	
基因 ID (Gene ID)	10246.0
蛋白 ID (SwissProt ID)	O00624.The antiserum was produced against synthesized peptide derived from human SLC17A2. AA range:213-262

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	47kDa

产品名称: SLC17A2 Rabbit Polyclonal Antibody
产品货号: APRab17953




研究背景 (Background)

function: May be involved in actively transporting phosphate into cells via Na(+) cotransport., similarity: Belongs to the major facilitator superfamily. Sodium/anion cotransporter family., tissue specificity: Expressed in the small intestine, kidney, spleen and testis. Not detected in fetal brain, bone marrow, and mammary gland., function: May be involved in actively transporting phosphate into cells via Na(+) cotransport., similarity: Belongs to the major facilitator superfamily. Sodium/anion cotransporter family., tissue specificity: Expressed in the small intestine, kidney, spleen and testis. Not detected in fetal brain, bone marrow, and mammary gland.,


研究领域 (Research Area)

Plasma Membrane; Channels; Signal Transduction; Metabolism; Vitamins / Minerals

图片 (Image Data)




Western blot analysis of lysates from HepG2 cells, using SLC17A2 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from RAW264.7 cells using SLC17A2 antibody.

产品名称: SLC17A2 Rabbit Polyclonal Antibody
产品货号: APRab17953



Western Blot analysis of various cells using SLC17A2 Polyclonal Antibody

注意事项 (Note)

For research use only .