

产品名称: Slit3 Rabbit Polyclonal Antibody
产品货号: AP Rab17974



产品概述 (Summary)

产品名称 (Production Name)	Slit3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SLIT3 SLIT3; KIAA0814; MEGF5; SLIL2; Slit homolog 3 protein; Slit-3; Multiple epidermal growth factor-like domains protein 5; Multiple EGF-like domains protein 5
别名 (Alternative Names)	
基因 ID (Gene ID)	6586.0
蛋白 ID (SwissProt ID)	O75094.Synthesized peptide derived from the Internal region of human Slit3.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	

产品名称: Slit3 Rabbit Polyclonal Antibody
产品货号: APRab17974




研究背景 (Background)

The protein encoded by this gene is secreted, likely interacting with roundabout homolog receptors to effect cell migration. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2012],function:May act as molecular guidance cue in cellular migration, and function may be mediated by interaction with roundabout homolog receptors.,sequence caution:Intron retention.,similarity:Contains 1 CTCK (C-terminal cystine knot-like) domain.,similarity:Contains 1 laminin G-like domain.,similarity:Contains 23 LRR (leucine-rich) repeats.,similarity:Contains 9 EGF-like domains.,tissue specificity:Predominantly expressed in thyroid.,

研究领域 (Research Area)

Axon guidance;

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200 (4° overnight) . 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200 (room temperature, 30min) .

注意事项 (Note)

For research use only .