

产品名称: SMF Rabbit Polyclonal Antibody
产品货号: APRab18018



产品概述 (Summary)

产品名称 (Production Name)	SMF Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	HMGXB3
别名 (Alternative Names)	HMGXB3; KIAA0194; SMF; HMG domain-containing protein 3; HMG box-containing protein 3; Protein SMF
基因 ID (Gene ID)	22993.0
蛋白 ID (SwissProt ID)	Q12766. The antiserum was produced against synthesized peptide derived from human HMGXB3. AA range: 131-180

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: SMF Rabbit Polyclonal Antibody
产品货号: APRab18018




研究背景 (Background)

This gene is one of the non-canonical high mobility group (HMG) genes. The encoded protein contains an HMG-box domain found in DNA binding proteins such as transcription factors and chromosomal proteins. [provided by RefSeq, Aug 2011],similarity:Contains 1 HMG box DNA-binding domain.,

研究领域 (Research Area)

Epigenetics and Nuclear Signaling; DNA / RNA; RNA Processing; Splicing

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma, using SMF Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .