

产品名称: SMVT Rabbit Polyclonal Antibody
产品货号: APRab18038



产品概述 (Summary)

产品名称 (Production Name)	SMVT Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SLC5A6
别名 (Alternative Names)	SLC5A6; SMVT; Sodium-dependent multivitamin transporter; Na(+)-dependent multivitamin transporter; Solute carrier family 5 member 6
基因 ID (Gene ID)	8884.0
蛋白 ID (SwissProt ID)	Q9Y289.The antiserum was produced against synthesized peptide derived from human SLC5A6. AA range:551-600

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	68kDa

产品名称: SMVT Rabbit Polyclonal Antibody
产品货号: APRab18038




研究背景 (Background)

function:Transports pantothenate, biotin and lipoate in the presence of sodium.,similarity:Belongs to the sodium:solute symporter (SSF) (TC 2.A.21) family.,function:Transports pantothenate, biotin and lipoate in the presence of sodium.,similarity:Belongs to the sodium:solute symporter (SSF) (TC 2.A.21) family.,


研究领域 (Research Area)

Plasma Membrane; Channels; Signal Transduction; Metabolism; Vitamins / Minerals; Pathways and Processes; Cofactors, Vitamins / minerals; Vitamins / minerals

图片 (Image Data)



Western blot analysis of lysates from HepG2 and HUVEC cells, using SLC5A6 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from COLO205 cells using SLC5A6 antibody.

注意事项 (Note)

For research use only .