

产品名称: SNX4 Rabbit Polyclonal Antibody
产品货号: APRab18077



产品概述 (Summary)

产品名称 (Production Name)	SNX4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SNX4
别名 (Alternative Names)	SNX4; Sorting nexin-4
基因 ID (Gene ID)	8723.0
蛋白 ID (SwissProt ID)	O95219.The antiserum was produced against synthesized peptide derived from human SNX4. AA range:191-240

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	60kDa

研究背景 (Background)

产品名称: SNX4 Rabbit Polyclonal Antibody
产品货号: APRab18077




This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein associated with the long isoform of the leptin receptor and with receptor tyrosine kinases for platelet-derived growth factor, insulin, and epidermal growth factor in cell cultures, but its function is unknown. This protein may form oligomeric complexes with family members. Two transcript variants, one protein coding and the other non-protein coding, have been found for this gene. [provided by RefSeq, Nov 2012],function:May be involved in several stages of intracellular trafficking.,similarity:Belongs to the sorting nexin family.,similarity:Contains 1 PX (phox homology) domain.,

研究领域 (Research Area)

Signal Transduction; Protein Trafficking; Vesicle Transport; Regulation

图片 (Image Data)



Western blot analysis of the lysates from HeLa cells using SNX4 antibody.

注意事项 (Note)

For research use only .