

产品名称: SPINK6 Rabbit Polyclonal Antibody
产品货号: APRab18197



产品概述 (Summary)

产品名称 (Production Name)	SPINK6 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SPINK6
别名 (Alternative Names)	SPINK6; Serine protease inhibitor Kazal-type 6; Kallikrein inhibitor
基因 ID (Gene ID)	404203.0
蛋白 ID (SwissProt ID)	Q6UWN8. The antiserum was produced against synthesized peptide derived from human SPINK6. AA range: 13-62

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)

产品名称: SPINK6 Rabbit Polyclonal Antibody
产品货号: APRab18197




The protein encoded by this gene is a Kazal-type serine protease inhibitor that acts on kallikrein-related peptidases in the skin. Two transcript variants the same protein have been found for this gene. [provided by RefSeq, Aug 2010],function:Probable serine protease inhibitor.,similarity:Contains 1 Kazal-like domain.,

研究领域 (Research Area)

Cell Biology; Proteolysis / Ubiquitin; Protease inhibitors; Other protease inhibitors

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma, using SPINK6 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .