

产品名称: SRPK3 Rabbit Polyclonal Antibody
产品货号: AP Rab18284



产品概述 (Summary)

产品名称 (Production Name)	SRPK3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SRPK3 SRPK3; MSK1; STK23; SRSF protein kinase 3; Muscle-specific serine kinase 1;
别名 (Alternative Names)	MSK-1; Serine/arginine-rich protein-specific kinase 3; SR-protein-specific kinase 3; Serine/threonine-protein kinase 23
基因 ID (Gene ID)	26576.0
蛋白 ID (SwissProt ID)	Q9UPE1.The antiserum was produced against synthesized peptide derived from human SRPK3. AA range:261-310

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	60kDa

产品名称: SRPK3 Rabbit Polyclonal Antibody
产品货号: APRab18284




研究背景 (Background)

This gene encodes a protein kinase similar to a protein kinase which is specific for the SR (serine/arginine-rich domain) family of splicing factors. A highly similar protein has been shown to play a role in muscle development in mice. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2009], catalytic activity: ATP + a protein = ADP + a phosphoprotein., function: Required for normal muscle development., similarity: Belongs to the protein kinase superfamily., similarity: Belongs to the protein kinase superfamily. CMGC Ser/Thr protein kinase family., similarity: Contains 1 protein kinase domain., tissue specificity: Exclusively expressed in skeletal and heart muscle.,


研究领域 (Research Area)

Signal Transduction; Protein Phosphorylation; Ser / Thr Kinases; Other Kinases; Developmental Biology; Organogenesis; Skeletal development; Muscle

图片 (Image Data)



Western blot analysis of SRPK3 Antibody. The lane on the right is blocked with the SRPK3 peptide.



Western blot analysis of the lysates from HeLa cells using SRPK3 antibody.

注意事项 (Note)

For research use only .

产品名称: SRPK3 Rabbit Polyclonal Antibody
产品货号: APRab18284

