

产品名称: STAC2 Rabbit Polyclonal Antibody
产品货号: AP Rab18330



产品概述 (Summary)

产品名称 (Production Name)	STAC2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	STAC2
别名 (Alternative Names)	STAC2; SH3 and cysteine-rich domain-containing protein 2; 24b2/STAC2; Src homology 3 and cysteine-rich domain-containing protein 2
基因 ID (Gene ID)	342667.0
蛋白 ID (SwissProt ID)	Q6ZMT1.The antiserum was produced against synthesized peptide derived from human STAC2. AA range:211-260

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	40kDa

产品名称: STAC2 Rabbit Polyclonal Antibody
产品货号: APRab18330




研究背景 (Background)

similarity:Contains 1 phorbol-ester/DAG-type zinc finger.,similarity:Contains 1 SH3 domain.,similarity:Contains 1 phorbol-ester/DAG-type zinc finger.,similarity:Contains 1 SH3 domain.,

研究领域 (Research Area)

Signal Transduction; Signaling Pathway; Nuclear Signaling; STATs; Epigenetics and Nuclear Signaling; Transcription; Microbiology; Organism; Virus; RNA Virus; ssRNA positive strand virus; SARS Coronavirus

图片 (Image Data)



Western blot analysis of lysates from K562 cells, using STAC2 Antibody. The lane on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .