

产品名称: STAC3 Rabbit Polyclonal Antibody
产品货号: APRab18331



产品概述 (Summary)

产品名称 (Production Name)	STAC3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	STAC3
别名 (Alternative Names)	
基因 ID (Gene ID)	246329.0
蛋白 ID (SwissProt ID)	Q96MF2.Synthesized peptide derived from part region of human protein

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	40kDa

研究背景 (Background)

The protein encoded by this gene is a component of the excitation-contraction coupling machinery of muscles. This protein is a member of the Stac gene family and contains an N-terminal cysteine-rich domain and two SH3 domains. Mutations in

产品名称: STAC3 Rabbit Polyclonal Antibody
产品货号: APRab18331




this gene are a cause of Native American myopathy. [provided by RefSeq, Nov 2013],similarity:Contains 1 phorbol-ester/DAG-type zinc finger.,similarity:Contains 2 SH3 domains.,

研究领域 (Research Area)

Cell Biology; Other Antibodies

图片 (Image Data)



Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4°over night

注意事项 (Note)

For research use only .