

产品名称: StARD10 Rabbit Polyclonal Antibody
产品货号: APRab18341



产品概述 (Summary)

产品名称 (Production Name)	StARD10 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	STARD10 STARD10; SDCCAG28; CGI-52; PCTP-like protein; PCTP-L; Antigen NY-CO-28;
别名 (Alternative Names)	START domain-containing protein 10; StARD10; Serologically defined colon cancer antigen 28; StAR-related lipid transfer protein 10
基因 ID (Gene ID)	10809.0
蛋白 ID (SwissProt ID)	Q9Y365.Synthesized peptide derived from the Internal region of human StARD10.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	33kDa

产品名称: StARD10 Rabbit Polyclonal Antibody
产品货号: APRab18341




研究背景 (Background)

function:May play specific roles in sperm maturation or fertilization.,sequence caution:Various sequencing problems as well as a translation in a wrong frame.,similarity:Contains 1 START domain.,function:May play specific roles in sperm maturation or fertilization.,sequence caution:Various sequencing problems as well as a translation in a wrong frame.,similarity:Contains 1 START domain.,

研究领域 (Research Area)

Signal Transduction; Metabolism; Lipid metabolism; Cancer; Tumor biomarkers

图片 (Image Data)



Western Blot analysis of CoLo cells using StARD10 Polyclonal Antibody

注意事项 (Note)

For research use only .