

产品名称: Striatin Rabbit Polyclonal Antibody
产品货号: AP Rab18408



产品概述 (Summary)

产品名称 (Production Name)	Striatin Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	STRN
别名 (Alternative Names)	STRN; Striatin
基因 ID (Gene ID)	6801.0
蛋白 ID (SwissProt ID)	O43815.Synthesized peptide derived from the C-terminal region of human Striatin.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	86kDa

研究背景 (Background)


产品名称: Striatin Rabbit Polyclonal Antibody
产品货号: AP Rab18408




function: Calmodulin-binding protein which may function as scaffolding or signaling protein and may play a role in dendritic Ca(2+) signaling., similarity: Belongs to the WD repeat striatin family., similarity: Contains 6 WD repeats., subunit: Interacts with protein phosphatase 2A (PP2A) ., tissue specificity: Preferentially expressed in brain., function: Calmodulin-binding protein which may function as scaffolding or signaling protein and may play a role in dendritic Ca(2+) signaling., similarity: Belongs to the WD repeat striatin family., similarity: Contains 6 WD repeats., subunit: Interacts with protein phosphatase 2A (PP2A) ., tissue specificity: Preferentially expressed in brain.,

研究领域 (Research Area)

图片 (Image Data)




Western Blot analysis of extracts from K562 cells, using Striatin Polyclonal Antibody.. Secondary antibody was diluted at 1:20000



Immunohistochemical analysis of paraffin-embedded mouse-brain, antibody was diluted at 1:100

产品名称: Striatin Rabbit Polyclonal Antibody
产品货号: APRab18408



Immunohistochemical analysis of paraffin-embedded mouse-brain, antibody was diluted at 1:100

注意事项 (Note)

For research use only .