

产品名称: STRN4 Rabbit Polyclonal Antibody
产品货号: APRab18410



产品概述 (Summary)

产品名称 (Production Name)	STRN4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	STRN4 ZIN
别名 (Alternative Names)	
基因 ID (Gene ID)	29888.0
蛋白 ID (SwissProt ID)	Q9NRL3.Synthesized peptide derived from human protein . at AA range: 250-330

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	82kDa

研究背景 (Background)

function: Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling

产品名称: STRN4 Rabbit Polyclonal Antibody
产品货号: APRab18410




protein.,miscellaneous:The name 'Zinedin' probably originates from the name of the famous soccer player from Marseille (Zinedine Zidane),similarity:Belongs to the WD repeat striatin family.,similarity:Contains 7 WD repeats.,function:Binds calmodulin in a calcium dependent manner. May function as scaffolding or signaling protein.,miscellaneous:The name 'Zinedin' probably originates from the name of the famous soccer player from Marseille (Zinedine Zidane),similarity:Belongs to the WD repeat striatin family.,similarity:Contains 7 WD repeats.,

研究领域 (Research Area)

Neuroscience

图片 (Image Data)



Western Blot analysis of Hela lysis, using primary antibody at 1:1000 dilution. Secondary antibody was diluted at 1:10000

注意事项 (Note)

For research use only .