

产品名称: SUHW4 Rabbit Polyclonal Antibody
产品货号: APRab18435



产品概述 (Summary)

产品名称 (Production Name)	SUHW4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZNF280D
别名 (Alternative Names)	ZNF280D; KIAA1584; SUHW4; ZNF634; Zinc finger protein 280D; Suppressor of hairy wing homolog 4; Zinc finger protein 634
基因 ID (Gene ID)	54816.0
蛋白 ID (SwissProt ID)	Q6N043.Synthesized peptide derived from the Internal region of human SUHW4.

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: SUHW4 Rabbit Polyclonal Antibody
产品货号: APRab18435




研究背景 (Background)

function:May function as a transcription factor.,similarity:Contains 5 C2H2-type zinc fingers.,function:May function as a transcription factor.,similarity:Contains 5 C2H2-type zinc fingers.,

研究领域 (Research Area)

Tags & Cell Markers; Subcellular Markers; Organelles; ER

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200 (4° overnight) . 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200 (room temperature, 45min) .

注意事项 (Note)

For research use only .