

产品名称: SVOP Rabbit Polyclonal Antibody
产品货号: APRab18467



产品概述 (Summary)

产品名称 (Production Name)	SVOP Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SVOP
别名 (Alternative Names)	SVOP; Synaptic vesicle 2-related protein; SV2-related protein
基因 ID (Gene ID)	55530.0
蛋白 ID (SwissProt ID)	Q8N4V2.The antiserum was produced against synthesized peptide derived from human SVOP. AA range:271-320

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	61kDa

研究背景 (Background)


产品名称: SVOP Rabbit Polyclonal Antibody
产品货号: APRab18467




similarity:Belongs to the major facilitator superfamily.,similarity:Belongs to the major facilitator superfamily.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from K562 cells, using SVOP Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from Jurkat cells using SVOP antibody.

注意事项 (Note)

For research use only .