

产品名称: SWAP Rabbit Polyclonal Antibody
产品货号: APRab18469



产品概述 (Summary)

产品名称 (Production Name)	SWAP Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SFSWAP SFSWAP; SFRS8; SWAP; Splicing factor; suppressor of white-apricot homolog;
别名 (Alternative Names)	Splicing factor, arginine/serine-rich 8; Suppressor of white apricot protein homolog
基因 ID (Gene ID)	6433.0
蛋白 ID (SwissProt ID)	Q12872.The antiserum was produced against synthesized peptide derived from human SFRS8. AA range:201-250

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	105kDa

产品名称: SWAP Rabbit Polyclonal Antibody
产品货号: APRab18469




研究背景 (Background)


This gene encodes a human homolog of Drosophila splicing regulatory protein. This gene autoregulates its expression by control of splicing of its first two introns. In addition, it also regulates the splicing of fibronectin and CD45 genes. Two transcript variants encoding different isoforms have been identified. [provided by RefSeq, May 2012],function:May function as an alternative splicing regulator. Regulate its own expression at the level of RNA processing. Also regulates the splicing of fibronectin and CD45 genes. May act, at least in part, by interaction with other R/S-containing splicing factors.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 2 SURP motif repeats.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of SFRS8 Antibody. The lane on the right is blocked with the SFRS8 peptide.



Western blot analysis of the lysates from 293 cells using SFRS8 antibody.

产品名称: SWAP Rabbit Polyclonal Antibody
产品货号: APRab18469



Western Blot analysis of various cells using SWAP Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .

注意事项 (Note)

For research use only .