

产品名称: T2R31 Rabbit Polyclonal Antibody
产品货号: APRab18564



产品概述 (Summary)

产品名称 (Production Name)	T2R31 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TAS2R31 TAS2R44
别名 (Alternative Names)	
基因 ID (Gene ID)	259290.0
蛋白 ID (SwissProt ID)	P59538.Synthesized peptide derived from human protein . at AA range: 10-90

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	33kDa

研究背景 (Background)

TAS2R44 belongs to the large TAS2R receptor family. TAS2Rs are expressed on the surface of taste receptor cells and mediate the perception of bitterness through a G protein-coupled second messenger pathway (Conte et al., 2002 [PubMed])

产品名称: T2R31 Rabbit Polyclonal Antibody
产品货号: AP Rab18564




12584440)). For further information on TAS2Rs, see MIM 604791.[supplied by OMIM, Mar 2009],cell surface receptor linked signal transduction, G-protein coupled receptor protein signaling pathway, sensory perception, sensory perception of chemical stimulus, neurological system process, cognition, sensory perception of taste,

研究领域 (Research Area)

Taste transduction;

图片 (Image Data)



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4°over night

注意事项 (Note)

For research use only .