

产品名称: T2R42 Rabbit Polyclonal Antibody
产品货号: APRab18569



产品概述 (Summary)

产品名称 (Production Name)	T2R42 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TAS2R42 TAS2R55
别名 (Alternative Names)	
基因 ID (Gene ID)	353164.0
蛋白 ID (SwissProt ID)	Q7RTR8.Synthesized peptide derived from human protein . at AA range: 110-190

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	34kDa

研究背景 (Background)

function:Receptor that may play a role in the perception of bitterness and is gustducin-linked. May play a role in sensing

产品名称: T2R42 Rabbit Polyclonal Antibody
产品货号: APRab18569




the chemical composition of the gastrointestinal content. The activity of this receptor may stimulate alpha gustducin, mediate PLC-beta-2 activation and lead to the gating of TRPM5.,miscellaneous:Most taste cells may be activated by a limited number of bitter compounds; individual taste cells can discriminate among bitter stimuli.,similarity:Belongs to the G-protein coupled receptor T2R family.,function:Receptor that may play a role in the perception of bitterness and is gustducin-linked. May play a role in sensing the chemical composition of the gastrointestinal content. The activity of this receptor may stimulate alpha gustducin, mediate PLC-beta-2 activation and lead to the gating of TRPM5.,miscellaneous:Most taste cells may be activated by a limited number of bitter compounds; individual taste cells can discriminate among bitter stimuli.,similarity:Belongs to the G-protein coupled receptor T2R family.,

研究领域 (Research Area)

Taste transduction;

图片 (Image Data)



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4°over night

注意事项 (Note)

For research use only .