

产品名称: T3JAM Rabbit Polyclonal Antibody
产品货号: APRab18581



产品概述 (Summary)

产品名称 (Production Name)	T3JAM Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TRAF3IP3
别名 (Alternative Names)	TRAF3IP3; T3JAM; TRAF3-interacting JNK-activating modulator; TRAF3-interacting protein 3
基因 ID (Gene ID)	80342.0
蛋白 ID (SwissProt ID)	Q9Y228.The antiserum was produced against synthesized peptide derived from human T3JAM. AA range:251-300

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	64kDa

产品名称: T3JAM Rabbit Polyclonal Antibody
产品货号: APRab18581




研究背景 (Background)


The gene encodes a protein that mediates cell growth by modulating the c-Jun N-terminal kinase signal transduction pathway. The encoded protein may also interact with a large multi-protein assembly containing the phosphatase 2A catalytic subunit. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2013],function:May function as an adapter molecule that regulates TRAF3-mediated JNK activation.,subunit: Binds to the isoleucine zipper of TRAF3 via its coiled-coil domain.,

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence analysis of HeLa cells, using T3JAM Antibody. The picture on the right is blocked with the synthesized



Western blot analysis of lysates from Jurkat cells, using T3JAM Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: T3JAM Rabbit Polyclonal Antibody
产品货号: APRab18581



Western Blot analysis of various cells using T3JAM Polyclonal Antibody

注意事项 (Note)

For research use only .