

产品名称: TA2R9 Rabbit Polyclonal Antibody
产品货号: APRab18586



产品概述 (Summary)

产品名称 (Production Name)	TA2R9 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TAS2R9
别名 (Alternative Names)	
基因 ID (Gene ID)	50835.0
蛋白 ID (SwissProt ID)	Q9NYW1.Synthesized peptide derived from human protein . at AA range: 90-170

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	34kDa

研究背景 (Background)

This gene product belongs to the family of candidate taste receptors that are members of the G-protein-coupled receptor

产品名称: TA2R9 Rabbit Polyclonal Antibody
产品货号: APRab18586




superfamily. These proteins are specifically expressed in the taste receptor cells of the tongue and palate epithelia. They are organized in the genome in clusters and are genetically linked to loci that influence bitter perception in mice and humans. In functional expression studies, they respond to bitter tastants. This gene maps to the taste receptor gene cluster on chromosome 12p13. [provided by RefSeq, Jul 2008],function:Gustducin-coupled receptor implicated in the perception of bitter compounds in the oral cavity and the gastrointestinal tract. Signals through PLCB2 and the calcium-regulated cation channel TRPM5.,miscellaneous:Several bitter taste receptors are expressed in a single taste receptor cell.,similarity:Belongs to the G-protein coupled receptor T2R family.,tissue specificity:Expressed in subsets of taste receptor cells of the tongue and palate epithelium and exclusively in gustducin-positive cells.,

研究领域 (Research Area)

Taste transduction;

图片 (Image Data)



Western blot analysis of lysates from PC12 cells, primary antibody was diluted at 1:1000, 4°over night

注意事项 (Note)

For research use only .