

产品名称: TAB182 Rabbit Polyclonal Antibody
产品货号: APRab18595



产品概述 (Summary)

产品名称 (Production Name)	TAB182 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TNKS1BP1
别名 (Alternative Names)	TNKS1BP1; KIAA1741; TAB182; 182 kDa tankyrase-1-binding protein
基因 ID (Gene ID)	85456.0
蛋白 ID (SwissProt ID)	Q9C0C2.The antiserum was produced against synthesized peptide derived from human TNKS1BP1. AA range:1601-1650

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	150kDa

研究背景 (Background)


产品名称: TAB182 Rabbit Polyclonal Antibody
产品货号: APRab18595




PTM:ADP-ribosylated by TNKS1 (in vitro),PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,subcellular location:Colocalizes with chromosomes during mitosis, and in the cytoplasm with cortical actin.,subunit:Binds to the ANK repeat domain of TNKS1 and TNKS2.,tissue specificity:Detected in testis, ovary, lung, skeletal muscle, heart, prostate and pancreas, and at very low levels in brain and peripheral blood leukocytes.,PTM:ADP-ribosylated by TNKS1 (in vitro),PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,subcellular location:Colocalizes with chromosomes during mitosis, and in the cytoplasm with cortical actin.,subunit:Binds to the ANK repeat domain of TNKS1 and TNKS2.,tissue specificity:Detected in testis, ovary, lung, skeletal muscle, heart, prostate and pancreas, and at very low levels in brain and peripheral blood leukocytes.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from 293 cells, using TNKS1BP1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using TAB182 Polyclonal Antibody diluted at 1: 1000

注意事项 (Note)

For research use only .

产品名称: TAB182 Rabbit Polyclonal Antibody
产品货号: APRab18595

