

产品名称: TAUT Rabbit Polyclonal Antibody
产品货号: APRab18660



产品概述 (Summary)

产品名称 (Production Name)	TAUT Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02%
储存溶液 (Buffer)	New type preservative N
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SLC6A6
别名 (Alternative Names)	SLC6A6; Sodium- and chloride-dependent taurine transporter; Solute carrier family 6 member 6
基因 ID (Gene ID)	6533.0
蛋白 ID (SwissProt ID)	P31641.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	Calculated MW: 70 kDa;

产品名称: TAUT Rabbit Polyclonal Antibody
产品货号: APRab18660




研究背景 (Background)

This gene encodes a multi-pass membrane protein that is a member of a family of sodium and chloride-ion dependent transporters. The encoded protein transports taurine and beta-alanine. There is a pseudogene for this gene on chromosome 21. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2013],

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using SLC6A6 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .