

产品名称: TBC1D3A/B/C Rabbit Polyclonal Antibody
产品货号: APRab18676



产品概述 (Summary)

产品名称 (Production Name)	TBC1D3A/B/C Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TBC1D3 TBC1D3; PRC17; TBC1D3A; TBC1 domain family member 3; Prostate cancer gene 17 protein; Protein TRE17-alpha; Rab GTPase-activating protein PRC17;
别名 (Alternative Names)	TBC1D3B; TBC1D3I; TBC1 domain family member 3B; TBC1D3C; TBC1D3D; TBC1 domain family member 3C
基因 ID (Gene ID)	414059/727735/84218/643947/414060/653380/654341/729873
蛋白 ID (SwissProt ID)	Q8IZP1/A6NDS4/Q6IPX1. The antiserum was produced against synthesized peptide derived from human TBC1D3. AA range: 481-530

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
----------------------	---

产品名称: TBC1D3A/B/C Rabbit Polyclonal Antibody
产品货号: APRab18676



蛋白分子量 (Molecular Weight) 62kDa


研究背景 (Background)

This gene represents one of a cluster of related genes found on chromosome 17. The proteins encoded by this gene family contain a TBC (Tre-2, Bub2p, and Cdc16p) domain and may be involved in GTPase signaling and vesicle trafficking.


[provided by RefSeq, Apr 2014],disease:May be involved in forms of prostate cancers.,function:Acts as a GTPase activating protein for RAB5. Does not act on RAB4 or RAB11.,similarity:Contains 1 Rab-GAP TBC domain.,tissue specificity:Expressed in liver, skeletal muscle, kidney, pancreas, spleen, testis, ovary, small intestine and peripheral blood leukocytes. Overexpressed in prostate cancers.,tissue specificity:Expressed in pancreas, thymus and testis.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from HeLa cells, using TBC1D3 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using TBC1D3A/B/C Polyclonal Antibody

产品名称: TBC1D3A/B/C Rabbit Polyclonal Antibody
产品货号: APRab18676



注意事项 (Note)

For research use only .