

产品名称: TBPL2 Rabbit Polyclonal Antibody
产品货号: APRab18695



产品概述 (Summary)

产品名称 (Production Name)	TBPL2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TBPL2
别名 (Alternative Names)	TBPL2; TBP2; TRF3; TATA box-binding protein-like protein 2; TBP-like protein 2; TATA box-binding protein-related factor 3; TBP-related factor 3
基因 ID (Gene ID)	387332.0
蛋白 ID (SwissProt ID)	Q6SJ96.The antiserum was produced against synthesized peptide derived from human TBPL2. AA range:81-130

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	41kDa

产品名称: TBPL2 Rabbit Polyclonal Antibody
产品货号: APRab18695




研究背景 (Background)

function: Transcription factor required in complex with TAF3 for the differentiation of myoblasts into myocytes. The complex replaces TFIID at specific promoters at an early stage in the differentiation process., similarity: Belongs to the TBP family., subcellular location: Present in the cytoplasm during cytokinesis., subunit: Interacts with TAF3., tissue specificity: Ubiquitously expressed in all tissues examined with highest levels in heart, lung, ovary, spleen and testes., function: Transcription factor required in complex with TAF3 for the differentiation of myoblasts into myocytes. The complex replaces TFIID at specific promoters at an early stage in the differentiation process., similarity: Belongs to the TBP family., subcellular location: Present in the cytoplasm during cytokinesis., subunit: Interacts with TAF3., tissue specificity: Ubiquitously expressed in all tissues examined with highest levels in heart, lung, ovary, spleen and testes.,


研究领域 (Research Area)

Basal transcription factors; Huntington's disease;

图片 (Image Data)




Western blot analysis of TBPL2 Antibody. The lane on the right is blocked with the TBPL2 peptide.




Immunohistochemistry analysis of paraffin-embedded human lung carcinoma, using TBPL2 Antibody. The lane on the right is blocked with the TBPL2 peptide.

产品名称: TBPL2 Rabbit Polyclonal Antibody
产品货号: APRab18695



Western blot analysis of the lysates from COLO205 cells using TBPL2 antibody.



Western Blot analysis of HeLa cells using TBPL2 Polyclonal Antibody

注意事项 (Note)

For research use only .