

产品名称: TCEAL3/5/6 Rabbit Polyclonal Antibody  
产品货号: APRab18724



## 产品概述 (Summary)

|                        |                                       |
|------------------------|---------------------------------------|
| 产品名称 (Production Name) | TCEAL3/5/6 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody            |
| 宿主 (Host)              | Rabbit                                |
| 应用 (Application)       | IHC, ICC, IF, ELISA                   |
| 种属反应性 (Reactivity)     | Human, Rat, Mouse                     |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                                                                                                                                                |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | TCEAL3/5/6<br>TCEAL3; MSTP072; Transcription elongation factor A protein-like 3; TCEA-like protein 3; Transcription elongation factor S-II protein-like 3; TCEAL5; Transcription elongation factor A protein-like 5; TCEA-like protein 5; Transcription elonga |
| 别名 (Alternative Names) |                                                                                                                                                                                                                                                                |
| 基因 ID (Gene ID)        | 85012.0                                                                                                                                                                                                                                                        |
| 蛋白 ID (SwissProt ID)   | Q969E4/Q5H9L2/Q6IPX3. The antiserum was produced against synthesized peptide derived from human TCEAL3. AA range:141-190                                                                                                                                       |

## 产品应用 (Application)

|                      |                                                           |
|----------------------|-----------------------------------------------------------|
| 稀释比 (Dilution Ratio) | IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:20000-1:40000 |
|----------------------|-----------------------------------------------------------|

**产品名称: TCEAL3/5/6 Rabbit Polyclonal Antibody**  
**产品货号: APRab18724**



---


**蛋白分子量 (Molecular Weight)**

### **研究背景 (Background)**

This gene encodes a member of the transcription elongation factor A (SII)-like (TCEAL) gene family. Members of this family contain TFA domains and may function as nuclear phosphoproteins that modulate transcription in a promoter context-dependent manner. Multiple family members are located on the X chromosome. Alternative splicing results in multiple transcript variants encoding a single isoform. [provided by RefSeq, Jul 2008],function:May be involved in transcriptional regulation.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the TFS-II family. TFA subfamily.,

### **研究领域 (Research Area)**

### **图片 (Image Data)**



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using TCEAL3/5/6 Antibody. The picture on the right is blocked with the synthesized peptide.

### **注意事项 (Note)**

For research use only .