

产品名称: TCEAL4 Rabbit Polyclonal Antibody
产品货号: APRab18725



产品概述 (Summary)

产品名称 (Production Name)	TCEAL4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TCEAL4
别名 (Alternative Names)	TCEAL4; NPD017; Transcription elongation factor A protein-like 4; TCEA-like protein 4; Transcription elongation factor S-II protein-like 4
基因 ID (Gene ID)	79921.0
蛋白 ID (SwissProt ID)	Q96EI5. The antiserum was produced against synthesized peptide derived from human TCEAL4. AA range: 166-215

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: TCEAL4 Rabbit Polyclonal Antibody
产品货号: APRab18725




研究背景 (Background)

This gene encodes a member of the transcription elongation factor A (SII)-like (TCEAL) gene family. This family is comprised of nuclear phosphoproteins that modulate transcription in a promoter context-dependent manner. Multiple family members are located on the X chromosome. Alternatively splicing results in multiple transcript variants. There is a pseudogene for this gene on chromosome 13. [provided by RefSeq, Apr 2015],function:May be involved in transcriptional regulation.,similarity:Belongs to the TFS-II family. TFA subfamily.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using TCEAL4 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .