

产品名称: TCEAL6 Rabbit Polyclonal Antibody
产品货号: APRab18726



产品概述 (Summary)

产品名称 (Production Name)	TCEAL6 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TCEAL6
别名 (Alternative Names)	TCEAL6; Transcription elongation factor A protein-like 6; TCEA-like protein 6; Transcription elongation factor S-II protein-like 6
基因 ID (Gene ID)	158931.0
蛋白 ID (SwissProt ID)	Q6IPX3.The antiserum was produced against synthesized peptide derived from human TCEAL6. AA range:61-110

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	25kDa

产品名称: TCEAL6 Rabbit Polyclonal Antibody
产品货号: AP Rab18726




研究背景 (Background)

function: May be involved in transcriptional regulation., PTM: Phosphorylated upon DNA damage, probably by ATM or ATR., similarity: Belongs to the TFS-II family. TFA subfamily., function: May be involved in transcriptional regulation., PTM: Phosphorylated upon DNA damage, probably by ATM or ATR., similarity: Belongs to the TFS-II family. TFA subfamily.,


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)




Western blot analysis of lysates from HUVEC, 293, HepG2, and Jurkat cells, using TCEAL6 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HUVEC cells using TCEAL6 antibody.

产品名称: TCEAL6 Rabbit Polyclonal Antibody
产品货号: APRab18726



Western Blot analysis of various cells using TCEAL6 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA) .

注意事项 (Note)

For research use only .