

产品名称: TCF-19 Rabbit Polyclonal Antibody
产品货号: APRab18730



产品概述 (Summary)

产品名称 (Production Name)	TCF-19 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TCF19
别名 (Alternative Names)	TCF19; SC1; Transcription factor 19; TCF-19; Transcription factor SC1
基因 ID (Gene ID)	6941.0
蛋白 ID (SwissProt ID)	Q9Y242.Synthesized peptide derived from the Internal region of human TCF-19.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	35kDa

研究背景 (Background)

产品名称: TCF-19 Rabbit Polyclonal Antibody
产品货号: APRab18730




This gene encodes a protein that contains a PHD-type zinc finger domain and likely functions as a transcription factor. The encoded protein plays a role proliferation and apoptosis of pancreatic beta cells. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016],developmental stage:Growth regulated.,function:Potential trans-activating factor that could play an important role in the transcription of genes required for the later stages of cell cycle progression.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 FHA domain.,similarity:Contains 1 PHD-type zinc finger.,

研究领域 (Research Area)

Epigenetics and Nuclear Signaling; Transcription; Other factors; Stem Cells; Signaling Pathways; Wnt Nuclear

图片 (Image Data)



Western Blot analysis of extracts from Jurkat cells, using TCF-19 Polyclonal Antibody.. Secondary antibody was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .

注意事项 (Note)

For research use only .