

产品名称: Tctex5 Rabbit Polyclonal Antibody  
产品货号: APRab18755



## 产品概述 (Summary)

|                        |                                   |
|------------------------|-----------------------------------|
| 产品名称 (Production Name) | Tctex5 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody        |
| 宿主 (Host)              | Rabbit                            |
| 应用 (Application)       | IHC, ICC, IF, ELISA               |
| 种属反应性 (Reactivity)     | Human, Mouse, Rat                 |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                                     |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | PPP1R11                                                                                                                                             |
| 别名 (Alternative Names) | PPP1R11; HCGV; TCTE5; Protein phosphatase 1 regulatory subunit 11; Hemochromatosis candidate gene V protein; HCG V; Protein phosphatase inhibitor 3 |
| 基因 ID (Gene ID)        | 6992.0                                                                                                                                              |
| 蛋白 ID (SwissProt ID)   | O60927. The antiserum was produced against synthesized peptide derived from human PPP1R11. AA range: 71-120                                         |

## 产品应用 (Application)

|                          |                                                           |
|--------------------------|-----------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:20000-1:40000 |
| 蛋白分子量 (Molecular Weight) |                                                           |

产品名称: Tctex5 Rabbit Polyclonal Antibody  
产品货号: APRab18755




## 研究背景 (Background)

This gene encodes a specific inhibitor of protein phosphatase-1 (PP1) with a differential sensitivity toward the metal-independent and metal-dependent forms of PP1. The gene is located within the major histocompatibility complex class I region on chromosome 6. [provided by RefSeq, Jul 2008],function:Inhibitor of protein phosphatase 1.,tissue specificity:Widely expressed.,

## 研究领域 (Research Area)

## 图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using PPP1R11 Antibody. The picture on the right is blocked with the synthesized peptide.

## 注意事项 (Note)

For research use only .