

产品名称: TDE1 Rabbit Polyclonal Antibody
产品货号: APRab18758



产品概述 (Summary)

产品名称 (Production Name)	TDE1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SERINC3
别名 (Alternative Names)	SERINC3; DIFF33; TDE1; SBB199; Serine incorporator 3; Tumor differentially expressed protein 1
基因 ID (Gene ID)	10955.0
蛋白 ID (SwissProt ID)	Q13530.The antiserum was produced against synthesized peptide derived from human SERC3. AA range:361-410

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	52kDa

产品名称: TDE1 Rabbit Polyclonal Antibody
产品货号: APRab18758




研究背景 (Background)

function: May be involved in cellular transformation., PTM: N-glycosylated., similarity: Belongs to the TDE1 family., tissue specificity: Ubiquitous. Expression levels were increased fourfold to tenfold in lung tumor tissues compared with normal pulmonary tissues., function: May be involved in cellular transformation., PTM: N-glycosylated., similarity: Belongs to the TDE1 family., tissue specificity: Ubiquitous. Expression levels were increased fourfold to tenfold in lung tumor tissues compared with normal pulmonary tissues.,


研究领域 (Research Area)

Cancer; Oncoproteins/suppressors; Oncoproteins; Tumor biomarkers

图片 (Image Data)




Immunofluorescence analysis of HepG2 cells, using SERC3 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from 293 cells, using SERC3 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: TDE1 Rabbit Polyclonal Antibody
产品货号: APRab18758



Western Blot analysis of various cells using TDE1 Polyclonal Antibody

注意事项 (Note)

For research use only .