

产品名称: TDE2L Rabbit Polyclonal Antibody
产品货号: APRab18759



产品概述 (Summary)

产品名称 (Production Name)	TDE2L Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SERINC2
别名 (Alternative Names)	SERINC2; TDE2L; FKSG84; Serine incorporator 2; Tumor differentially expressed protein 2-like
基因 ID (Gene ID)	347735.0
蛋白 ID (SwissProt ID)	Q96SA4.The antiserum was produced against synthesized peptide derived from human SERC2. AA range:341-390

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	55kDa

产品名称: TDE2L Rabbit Polyclonal Antibody
产品货号: APRab18759




研究背景 (Background)


similarity:Belongs to the TDE1 family.,similarity:Belongs to the TDE1 family.,

研究领域 (Research Area)


图片 (Image Data)



Immunofluorescence analysis of NIH/3T3 cells, using SERC2 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from HeLa cells, using SERC2 Antibody. The lane on the right is blocked with the synthesized peptide.



产品名称: TDE2L Rabbit Polyclonal Antibody
产品货号: APRab18759



Western Blot analysis of various cells using TDE2L Polyclonal Antibody

注意事项 (Note)

For research use only .