

产品名称: TEF-4 Rabbit Polyclonal Antibody
产品货号: APRab18775



产品概述 (Summary)

产品名称 (Production Name)	TEF-4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TEAD2
别名 (Alternative Names)	TEAD2; TEF4; Transcriptional enhancer factor TEF-4; TEA domain family member 2; TEAD-2
基因 ID (Gene ID)	8463.0
蛋白 ID (SwissProt ID)	Q15562.The antiserum was produced against synthesized peptide derived from human TEAD2. AA range:71-120

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: TEF-4 Rabbit Polyclonal Antibody
产品货号: APRab18775




研究背景 (Background)


function: Putative transcription factor that binds to the SPH and GT-IIC "enhancers" (5'-GTGGAATGT-3'). May be involved in the gene regulation of neural development. Binds to the M-CAT motif., similarity: Contains 1 TEA DNA-binding domain., function: Putative transcription factor that binds to the SPH and GT-IIC "enhancers" (5'-GTGGAATGT-3'). May be involved in the gene regulation of neural development. Binds to the M-CAT motif., similarity: Contains 1 TEA DNA-binding domain.,

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of HepG2 cells, using TEAD2 Antibody. The picture on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using TEAD2 Antibody. The picture on the right is blocked with the synthesized peptide.

产品名称: TEF-4 Rabbit Polyclonal Antibody
产品货号: APRab18775



Western Blot analysis of mouse-brain cells using TEF-4 Polyclonal Antibody cells nucleus extracted by Minute TM
Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .

注意事项 (Note)

For research use only .