

产品名称: TEF-5 Rabbit Polyclonal Antibody
产品货号: APRab18776



产品概述 (Summary)

产品名称 (Production Name)	TEF-5 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TEAD3
别名 (Alternative Names)	TEAD3; TEAD5; TEF5; Transcriptional enhancer factor TEF-5; DTEF-1; TEA domain family member 3; TEAD-3
基因 ID (Gene ID)	7005.0
蛋白 ID (SwissProt ID)	Q99594.The antiserum was produced against synthesized peptide derived from human TEAD3. AA range:211-260

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	49kDa

产品名称: TEF-5 Rabbit Polyclonal Antibody
产品货号: APRab18776




研究背景 (Background)

This gene product is a member of the transcriptional enhancer factor (TEF) family of transcription factors, which contain the TEA/ATTS DNA-binding domain. It is predominantly expressed in the placenta and is involved in the transactivation of the chorionic somatomammotropin-B gene enhancer. Translation of this protein is initiated at a non-AUG (AUA) start codon. [provided by RefSeq, Jul 2008],function: Binds to multiple functional elements of the human chorionic somatomammotropin-B gene enhancer.,similarity: Contains 1 TEA DNA-binding domain.,tissue specificity: Preferentially expressed in the placenta.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of the lysates from HUVEC cells using TEAD3 antibody.

注意事项 (Note)

For research use only .