

产品名称: Tetranectin Rabbit Polyclonal Antibody
产品货号: AP Rab18806



产品概述 (Summary)

产品名称 (Production Name)	Tetranectin Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CLEC3B TNA
别名 (Alternative Names)	Tetranectin (TN; C-type lectin domain family 3 member B; Plasminogen kringle 4-binding protein)
基因 ID (Gene ID)	7123.0
蛋白 ID (SwissProt ID)	P05452. Synthetic peptide from human protein at AA range: 11-60

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:200, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)

产品名称: Tetranectin Rabbit Polyclonal Antibody
产品货号: APRab18806




function:Tetranectin binds to plasminogen and to isolated kringle 4. May be involved in the packaging of molecules destined for exocytosis.,mass spectrometry: PubMed:10614823,online information:Tetranectin,similarity:Contains 1 C-type lectin domain.,subunit:Homotrimer.,tissue specificity:Found in plasma.,function:Tetranectin binds to plasminogen and to isolated kringle 4. May be involved in the packaging of molecules destined for exocytosis.,mass spectrometry: PubMed:10614823,online information:Tetranectin,similarity:Contains 1 C-type lectin domain.,subunit:Homotrimer.,tissue specificity:Found in plasma.,


研究领域 (Research Area)

Cardiovascular

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200

注意事项 (Note)

For research use only .