

产品名称: TFIIH p62 Rabbit Polyclonal Antibody
产品货号: APRab18833



产品概述 (Summary)

产品名称 (Production Name)	TFIIH p62 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	GTF2H1 GTF2H1; BTF2; General transcription factor IIH subunit 1; Basic transcription factor 2 62 kDa subunit; BTF2 p62; General transcription factor IIH polypeptide 1; TFIIH basal transcription factor complex p62 subunit
别名 (Alternative Names)	
基因 ID (Gene ID)	2965.0
蛋白 ID (SwissProt ID)	P32780.The antiserum was produced against synthesized peptide derived from human TF2H1. AA range:15-64

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	57kDa

产品名称: TFIIH p62 Rabbit Polyclonal Antibody
产品货号: AP Rab18833


研究背景 (Background)

function: Component of the core-TFIIH basal transcription factor involved in nucleotide excision repair (NER) of DNA and, when complexed to CAK, in RNA transcription by RNA polymerase II., PTM: Phosphorylated., similarity: Contains 2 BSD domains., subunit: One of the six subunits forming the core-TFIIH basal transcription factor. Interacts with PUF60., function: Component of the core-TFIIH basal transcription factor involved in nucleotide excision repair (NER) of DNA and, when complexed to CAK, in RNA transcription by RNA polymerase II., PTM: Phosphorylated., similarity: Contains 2 BSD domains., subunit: One of the six subunits forming the core-TFIIH basal transcription factor. Interacts with PUF60.,


研究领域 (Research Area)

Basal transcription factors; Nucleotide excision repair;

图片 (Image Data)




Immunohistochemistry analysis of paraffin-embedded human brain tissue, using TF2H1 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from Jurkat cells, using TF2H1 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: TFIIH p62 Rabbit Polyclonal Antibody
产品货号: APRab18833



Western blot analysis of the lysates from HepG2 cells using TF2H1 antibody.

注意事项 (Note)

For research use only .