

产品名称: THAP4 Rabbit Polyclonal Antibody
产品货号: AP Rab18872



产品概述 (Summary)

产品名称 (Production Name)	THAP4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	THAP4
别名 (Alternative Names)	THAP4; CGI-36; PP238; THAP domain-containing protein 4
基因 ID (Gene ID)	51078.0
蛋白 ID (SwissProt ID)	Q8WY91.The antiserum was produced against synthesized peptide derived from human THAP4. AA range:341-390

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	63kDa

研究背景 (Background)


产品名称: THAP4 Rabbit Polyclonal Antibody
产品货号: AP Rab18872




PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 THAP-type zinc finger.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 THAP-type zinc finger.,

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of lysates from COS7 cells, using THAP4 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HT-29 cells using THAP4 antibody.



Western Blot analysis of various cells using THAP4 Polyclonal Antibody

产品名称: THAP4 Rabbit Polyclonal Antibody
产品货号: APRab18872



注意事项 (Note)

For research use only .