

产品名称: TIF-IA Rabbit Polyclonal Antibody
产品货号: AP Rab18930



产品概述 (Summary)

产品名称 (Production Name)	TIF-IA Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RRN3
别名 (Alternative Names)	RNA polymerase I-specific transcription initiation factor RRN3 (Transcription initiation factor IA) (TIF-IA)
基因 ID (Gene ID)	54700.0
蛋白 ID (SwissProt ID)	Q9NYV6.Synthesized peptide derived from human TIF-IA Polyclonal

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	70kDa

研究背景 (Background)

产品名称: TIF-IA Rabbit Polyclonal Antibody
产品货号: APRab18930




function:Required for efficient transcription initiation by RNA polymerase I.,similarity:Belongs to the RRN3 family.,function:Required for efficient transcription initiation by RNA polymerase I.,similarity:Belongs to the RRN3 family.,

研究领域 (Research Area)

Epigenetics and Nuclear Signaling; Transcription; Polymerase associated factors; Pol I Transcription; Other factors

图片 (Image Data)



Western blot analysis of 3T3 lysate, antibody was diluted at 1000. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .