

产品名称: TIS11D Rabbit Polyclonal Antibody
产品货号: APRab18968



产品概述 (Summary)

产品名称 (Production Name)	TIS11D Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZFP36L2 ZFP36L2; BRF2; ERF2; RNF162C; TIS11D; Zinc finger protein 36; C3H1 type-like
别名 (Alternative Names)	2; ZFP36-like 2; Butyrate response factor 2; EGF-response factor 2; ERF-2; Protein TIS11D
基因 ID (Gene ID)	678.0
蛋白 ID (SwissProt ID)	P47974.The antiserum was produced against synthesized peptide derived from human TISD. AA range:207-256

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	55kDa

产品名称: TIS11D Rabbit Polyclonal Antibody
产品货号: APRab18968




研究背景 (Background)

This gene is a member of the TIS11 family of early response genes. Family members are induced by various agonists such as the phorbol ester TPA and the polypeptide mitogen EGF. The encoded protein contains a distinguishing putative zinc finger domain with a repeating cys-his motif. This putative nuclear transcription factor most likely functions in regulating the response to growth factors. [provided by RefSeq, Jul 2008],function:Probable regulatory protein involved in regulating the response to growth factors. RNA-binding protein that binds to 5'UUAUUUAUUU-3' core sequence. Binds to the class II AU-rich element (ARE) in the 3'-UTR of target mRNAs and promotes their deadenylation and degradation.,similarity:Contains 2 C3H1-type zinc fingers.,


研究领域 (Research Area)

Cell Biology; Apoptosis; Intracellular; p53 Pathway; Epigenetics and Nuclear Signaling; Transcription; Domain Families; Zinc Finger; Chromatin Binding Proteins; DNA / RNA binding

图片 (Image Data)



Western blot analysis of lysates from A549 cells, using TISD Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using TIS11D Polyclonal Antibody. Secondary antibody was diluted at 1:20000

产品名称: TIS11D Rabbit Polyclonal Antibody
产品货号: APRab18968



注意事项 (Note)

For research use only .