

产品名称: TM4SF19 Rabbit Polyclonal Antibody
产品货号: APRab19028



产品概述 (Summary)

产品名称 (Production Name)	TM4SF19 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TM4SF19 OCTM4
别名 (Alternative Names)	Transmembrane 4 L6 family member 19 (Osteoclast maturation-associated gene 4 protein;Tetraspan membrane protein OCTM4)
基因 ID (Gene ID)	116211.0
蛋白 ID (SwissProt ID)	Q96DZ7.Synthesized peptide derived from human TM4SF19. AA range:140-160

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	24kDa

产品名称: TM4SF19 Rabbit Polyclonal Antibody
产品货号: APRab19028




研究背景 (Background)

similarity:Belongs to the L6 tetraspanin family.,similarity:Belongs to the L6 tetraspanin family.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of various lysates, primary antibody was diluted at 1:1000, 4° over night, secondary antibody was diluted at 1:10000, 37° 1hour.

注意事项 (Note)

For research use only .