

产品名称: TMEM185A Rabbit Polyclonal Antibody
产品货号: APRab19053



产品概述 (Summary)

产品名称 (Production Name)	TMEM185A Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TMEM185A
别名 (Alternative Names)	TMEM185A; CXorf13; FAM11A; Transmembrane protein 185A; Protein FAM11A
基因 ID (Gene ID)	84548.0
蛋白 ID (SwissProt ID)	Q8NFB2.The antiserum was produced against synthesized peptide derived from human TMEM185A. AA range:290-339

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	41kDa

产品名称: TMEM185A Rabbit Polyclonal Antibody
产品货号: APRab19053




研究背景 (Background)


The protein encoded by this gene is predicted to be a transmembrane protein. This gene is best known for localizing to the CpG island of the fragile site FRAXF. The 5' untranslated region of this gene contains a CGG trinucleotide repeat sequence that normally consists of 7-40 tandem CGG repeats but which can expand to greater than 300 repeats. Methylation of the CpG island leads to transcriptional silencing of this gene, but neither the silencing nor an expanded repeat region appear to manifest itself in a clear phenotypic manner. Alternative splicing results in multiple transcript variants. A pseudogene of this gene has been defined on the X chromosome. [provided by RefSeq, Aug 2013],similarity:Belongs to the TMEM185 family.,

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of MCF7 cells, using TMEM185A Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from K562 cells, using TMEM185A Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: TMEM185A Rabbit Polyclonal Antibody
产品货号: APRab19053



Western blot analysis of the lysates from HT-29 cells using TMEM185A antibody.

注意事项 (Note)

For research use only .