

产品名称: Tmem72 Rabbit Polyclonal Antibody
产品货号: APRab19058



产品概述 (Summary)

产品名称 (Production Name)	Tmem72 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TMEM72 C10orf127 KSP37
别名 (Alternative Names)	Transmembrane protein 72 (Kidney-specific secretory protein of 37 kDa)
基因 ID (Gene ID)	643236.0
蛋白 ID (SwissProt ID)	A0PK05.Synthesized peptide derived from human Tmem72 Polyclonal

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	30kDa

研究背景 (Background)

This gene encodes a transmembrane protein which may be expressed specifically in the kidney. [provided by RefSeq, Sep

产品名称: Tmem72 Rabbit Polyclonal Antibody
产品货号: APRab19058




2016],

研究领域 (Research Area)

Cell Biology; Other Antibodies

图片 (Image Data)



Western blot analysis of various lysate, antibody was diluted at 1000. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .