

产品名称: TMPO Rabbit Polyclonal Antibody
产品货号: APRab19070



产品概述 (Summary)

产品名称 (Production Name)	TMPO Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TMPO LAP2
别名 (Alternative Names)	TMPO LAP2
基因 ID (Gene ID)	7112.0
蛋白 ID (SwissProt ID)	P42166/P42167.Synthetic peptide from human protein at AA range: 1-50

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	75+50kDa

研究背景 (Background)

The protein encoded by TMPO (thymopoietin) resides in the nucleus and may play a role in the assembly of the nuclear

产品名称: TMPO Rabbit Polyclonal Antibody
产品货号: APRab19070




lamina, and thus help maintain the structural organization of the nuclear envelope. It may function as a receptor for the attachment of lamin filaments to the inner nuclear membrane. Mutations in TMPO are associated with dilated cardiomyopathy. Alternatively spliced transcript variants encoding different isoforms have been noted for TMPO. regulation of transcription,


研究领域 (Research Area)

Signal Transduction; Cytoskeleton / ECM; Cytoskeleton; Intermediate Filaments; Class V; Lamins; Basal Lamina

图片 (Image Data)




Western Blot analysis of mouse-kidney, hepg2 cells using Antibody diluted at 800. Secondary antibody was diluted at 1:20000




Immunohistochemical analysis of paraffin-embedded human-cervical-cancer, antibody was diluted at 1:200


产品名称: TMPO Rabbit Polyclonal Antibody
产品货号: APRab19070



Immunohistochemical analysis of paraffin-embedded human-cervical-cancer, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200

注意事项 (Note)

For research use only .