

产品名称: TNF- β Rabbit Polyclonal Antibody
产品货号: APRab19097



产品概述 (Summary)

产品名称 (Production Name)	TNF- β Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LTA
别名 (Alternative Names)	LTA; TNFB; TNFSF1; Lymphotoxin-alpha; LT-alpha; TNF-beta; Tumor necrosis factor ligand superfamily member 1
基因 ID (Gene ID)	4049.0
蛋白 ID (SwissProt ID)	P01374.The antiserum was produced against synthesized peptide derived from human TNF-beta. AA range:17-66

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	22kDa

产品名称: TNF- β Rabbit Polyclonal Antibody
产品货号: APRab19097




研究背景 (Background)

The encoded protein, a member of the tumor necrosis factor family, is a cytokine produced by lymphocytes. The protein is highly inducible, secreted, and forms heterotrimers with lymphotoxin-beta which anchor lymphotoxin-alpha to the cell surface. This protein also mediates a large variety of inflammatory, immunostimulatory, and antiviral responses, is involved in the formation of secondary lymphoid organs during development and plays a role in apoptosis. Genetic variations in this gene are associated with susceptibility to leprosy type 4, myocardial infarction, non-Hodgkin's lymphoma, and psoriatic arthritis. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Jul 2012],disease:Genetic variations in LTA are associated with susceptibility to leprosy type 4 (LPRS4) [MIM:610988]; also known as susceptibility to age-dependent leprosy. Leprosy is a chronic infectious disease of peripheral sensory nerves caused by Mycobacterium leprae.,disease:Genetic variations in LTA are associated with susceptibility to psoriatic arthritis [MIM:607507]. Psoriasis is a chronic inflammatory dermatosis that affects approximately 2% of the population. It is characterized by red, scaly skin lesions that are usually found on the scalp, elbows, and knees, and may be associated with severe arthritis. Psoriatic arthritis has been defined as an inflammatory arthritis usually without any rheumatoid factor in serum (seronegative arthritis) associated with psoriasis.,function:Cytokine that in its homotrimeric form binds to TNFRSF1A/TNFR1, TNFRSF1B/TNFR and TNFRSF14/HVEM. In its heterotrimeric form with LTB binds to TNFRSF3/LTBR. Lymphotoxin is produced by lymphocytes and cytotoxic for a wide range of tumor cells in vitro and in vivo.,online information:Lymphotoxin entry,online information:The Singapore human mutation and polymorphism database,similarity:Belongs to the tumor necrosis factor family.,subcellular location:The homotrimer is secreted. The heterotrimer is membrane-associated.,subunit:Homotrimer, and heterotrimer of either two LTB and one LTA subunits or (less prevalent) two LTA and one LTB subunits.,

研究领域 (Research Area)


Cytokine-cytokine receptor interaction;Antigen processing and presentation;Type I diabetes mellitus;

图片 (Image Data)



Western blot analysis of lysate from COS-7 cell, using TNF- β antibody.

产品名称: TNF- β Rabbit Polyclonal Antibody
产品货号: APRab19097



Western Blot analysis of various cells using TNF- β Polyclonal Antibody diluted at 1: 1000. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .