

产品名称: TP1 Rabbit Polyclonal Antibody
产品货号: APRab19138



产品概述 (Summary)

产品名称 (Production Name)	TP1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TEP1
别名 (Alternative Names)	TEP1; TLP1; TP1; Telomerase protein component 1; Telomerase-associated protein 1; Telomerase protein 1; p240; p80 telomerase homolog
基因 ID (Gene ID)	7011.0
蛋白 ID (SwissProt ID)	Q99973. The antiserum was produced against synthesized peptide derived from human TEP1. AA range: 171-220

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: TP1 Rabbit Polyclonal Antibody
产品货号: APRab19138




研究背景 (Background)

This gene product is a component of the ribonucleoprotein complex responsible for telomerase activity which catalyzes the addition of new telomeres on the chromosome ends. The telomerase-associated proteins are conserved from ciliates to humans. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016],function:Component of the telomerase ribonucleoprotein complex that is essential for the replication of chromosome termini. Also component of the ribonucleoprotein vaults particle, a multi-subunit structure involved in nucleo-cytoplasmic transport. Responsible for the localizing and stabilizing vault RNA (vRNA) association in the vault ribonucleoprotein particle. Binds to TERC.,similarity:Contains 1 NACHT domain.,similarity:Contains 1 TROVE domain.,similarity:Contains 21 WD repeats.,similarity:Contains 4 TEP1 N-terminal repeats.,subunit:Component of the telomerase holoenzyme complex at least composed of TERT, DKC1, WDR79/TCAB1, NOP10, NHP2, GAR1, TEP1, EST1A, POT1 and a telomerase RNA template component (TERC). Component of the vault ribonucleoprotein particle, at least composed of MVP, PARP4 and one or more vault RNAs (vRNAs). Binds to VAULTRC1, VAULTRC2 and VAULTRC4/hvg4 vRNAs.,tissue specificity:Ubiquitous.,

研究领域 (Research Area)

Epigenetics and Nuclear Signaling; Transcription; Transcription Factors; Developmental Biology; Reproduction; Fertilisation; Histones; Variants

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma, using TEP1 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .